

BD Symphony A1

Laser	Channel	LP	BP	Fluorochromes
Violet (405 nm)	A	750	780/60	<ul style="list-style-type: none"> BV786 Qdot 800
	B	690	710/50	<ul style="list-style-type: none"> BV711 Qdot 700
	C	655	670/30	<ul style="list-style-type: none"> BV650 Qdot 655
	D	595	610/20	<ul style="list-style-type: none"> BV605 Qdot 605
	E	505	525/50	<ul style="list-style-type: none"> AmCyan BV510 V500 CFP Zombie Aqua
	F	410 410	431/28 450/50	<ul style="list-style-type: none"> PB450 BV421 Pacific Blue eFluor 450 DAPI
Blue (488 nm)	A	690	710/50	<ul style="list-style-type: none"> PerCP/Cy5.5 7-AAD
	B	505	530/30	<ul style="list-style-type: none"> FITC GFP Alexa488 CFSE MitoTracker Green
	C		488/10	<ul style="list-style-type: none"> SSC
Yellow Green (561 nm)	A	750	780/60	<ul style="list-style-type: none"> PE/Cy7
	B	690	710/50	<ul style="list-style-type: none"> PE/Cy5.5
	C	595	610/20	<ul style="list-style-type: none"> PE-TexasRed PI mCherry Alexa568 MitoTracker Red
	D	580	586/15	<ul style="list-style-type: none"> PE Alexa555 DsRed tdTomato mRFP MitoTracker Orange
Red (640 nm)	A	750	780/60	<ul style="list-style-type: none"> APC/Cy7 APC/efluor780 Alexa780 Live-Dead Near IR
	B	655	670/30	<ul style="list-style-type: none"> APC Alexa647 DRAQ5 / DRAQ7 Cy5 MitoTracker DeepRed