

BD ARIA Fusion

Laser	Channel	LP	BP	Fluorochromes
Violet (405 nm)	A	750	780/60	<ul style="list-style-type: none"> BV786 Qdot 800
	B	690	710/50	<ul style="list-style-type: none"> BV711 Qdot 700
	C	630	660/20	<ul style="list-style-type: none"> BV650 Qdot 655
	D	595	610/20	<ul style="list-style-type: none"> BV605 Qdot 605
	E	505	525/50	<ul style="list-style-type: none"> AmCyan BV510 V500 CFP Zombie Aqua
	F		450/50	<ul style="list-style-type: none"> PB450 BV421 Pacific Blue eFluor 450 DAPI
Blue (488 nm)	A	655	695/40	<ul style="list-style-type: none"> PerCP PerCP/Cy5.5
	B	502	530/30	<ul style="list-style-type: none"> FITC GFP Alexa488 CFSE MitoTracker Green
	C		488/10	<ul style="list-style-type: none"> SSC
Yellow Green (561 nm)	A	735	780/60	<ul style="list-style-type: none"> PE/Cy7
	B	685	710/50	<ul style="list-style-type: none"> PE/Cy5.5
	C	630	670/14	<ul style="list-style-type: none"> PE-Cy5 7-AAD MitoTracker Red FM
	D	600	610/20	<ul style="list-style-type: none"> PE-TexasRed PI mCherry Alexa568 MitoTracker Red
	E		582/15	<ul style="list-style-type: none"> PE Alexa555 DsRed tdTomato mRFP MitoTracker Orange
Red (640 nm)	A	755	780/60	<ul style="list-style-type: none"> APC/Cy5.5 APC/Cy7 APC/eFluor780 Alexa780 Live-Dead Near IR
	B	690	730/45	<ul style="list-style-type: none"> Alexa 700 DRAQ5 / DRAQ7
	C		670/30	<ul style="list-style-type: none"> APC Alexa647 DRAQ5 / DRAQ7 Cy5 MitoTracker DeepRed

