

CytoFLEX

Laser	Channel	BP	Fluorochromes
Violet (405 nm)	PB450	450/45	<ul style="list-style-type: none"> • PB450 • BV421 • Pacific Blue • eFluor 450 • DAPI
	KO525	525/40	<ul style="list-style-type: none"> • Krome Orange • BV510 • V500 • CFP • Zombie Aqua
	Violet610	610/20	<ul style="list-style-type: none"> • BV605 • Qdot 605
	Violet660	660/20	<ul style="list-style-type: none"> • BV650 • Qdot 655
	Violet780	780/60	<ul style="list-style-type: none"> • BV786 • Qdot 800
Blue (488 nm)	FITC	525/40	<ul style="list-style-type: none"> • FITC • Brilliant Blue 515 • Alexa488 • GFP • CFSE • Fluo-3
	PE	585/42	<ul style="list-style-type: none"> • PE • PI
	ECD	610/20	<ul style="list-style-type: none"> • ECD • PE-Texas Red • PE-CF594 • PE-Dazzle 594 • PI
	PC5.5	690/50	<ul style="list-style-type: none"> • PerCP, PerCP/Cy5.5 • PE/Cy5.5 • PI • 7-AAD
	PC7	780/60	<ul style="list-style-type: none"> • PE/Cy7 • PE/Alexa750
Red (633 nm)	APC	660/20	<ul style="list-style-type: none"> • APC • Alexa647 • eFluor 660
	APC/A700	712/25	<ul style="list-style-type: none"> • APC/A700 • BD Horizon APC/R700 • Alexa700
	APC/A750	780/60	<ul style="list-style-type: none"> • APC/A750 • APC-H7 • APC/Cy7 • APC/eFluor 780 • DRAQ5 • Zombie NIR